

WATER PROVIDES US LIFE



Are you sure that the water you drink is safe?



The drinking water is guaranteed as the output of the drinking water treatment plant and into the public taps, but not in private taps as your home, offices, factories, restaurants, hotels, schools, nursing homes, wells, reservoirs...



Public works in the street seriously endanger the security of the water.

Verify that there is no infiltration of sewage or natural water in your building network.



UNIQUE IN THE MARKET! Now you can check yourself with

DRINK WATER KIT

Complete kit for direct analysis of bacterial counts and presence / absence of indicators for wastewater infiltration: *Coliforms-E.coli*, faecal enterococci and *Clostridium perfringens* and its spores (European Directive on human use water).



In only 24-72 hours you will have the most reliable results by simple color changes of the 5 components of the kit.

- ★ Reference: KMT003, box for 10 complete test of 5 parameters each, in sampler bottles, for field or laboratory analysis
- ★ Also available in miniportable boxes for field analysis: MINI-DRINK WATER KIT Reference: KMT005, box for 3 complete test of 4 parameters each, in sampler bags

Tel: +34 918974616 - Fax: +34 918974641
P.O. Box 44, 28210 Madrid - Spain
Email: microkit@microkit.es www.microkit.es

