



5.

MICROKIT STAR PRODUCTS BY SECTORS

FOOD INDUSTRY AND CATERING, SUPERMARKETS:



* KIT PRO for immediate hygiene surfaces and water control KP030, 36 u

* CRIOTECA MIX to stock strains for years KPX200, 50u

* M-Ident Staph latex coagulase confirmation in colonies KWD094, 100u



* M-Ident Salmonella direct screening látex in enrichment broth KMB501, 50 u




* M-Ident *Clostridium perfringens* confirmation kit ISO 6461, KMT008, 9 test

* More than 200 references in dehydrated media, 100-500 grs
* More than 200 references in ready to use bottles & tubes



* ANAEROKIT, generator of anaerobic atmosphere, great bag for more plates. KKT001, 25 u



* Cromogenic PCA to diferenciate red colonies from sample particles, BCD510, 500 g for 1500 plates

* Bottles 200 ml Cromogenic PCA-Milk (red colonies), RPL215, 25 or 75 u (for 1000 plates) or dehydrated 100-500 grs



* Orange Serum Agar, total plate count in juices/marmalades, DMT092, 500 g

* Tributyrin Agar for lypolitic flora detection DMT125 and supplement SDA078



* MUGPLUS® MICROKIT® CROMOGENIC AGAR (BASE), DMT400 Coliform (Pink) and *E.coli* (Blue) in water (M.F.) and food 100-500 g + SMS400 CEFSULODIN/VANCOMYCIN sterile selective supplement 100 ml



* CHROMOSALM MICROKIT® AGAR Selective and differential isolement of Salmonella without interferences of Proteus, Citrobacter... DMT500, 100-500 g.

* CHROMOCYTOGENES MICROKIT® AGAR (BASE), Selective and differential isolement of *Listeria monocytogenes* ISO 11290 (15-10-2004), DMT700, 100-500 g, SMT700 Double Supplement Selective for Chromocytogenes Agar c.s.p.10 x 500 ml.



* M-IDENT® Rhamnose/Xylose *L.monocytogenes*, confirmation of green-with-halo colonies in Chromocytogenes Agar, KMT012, 10 or 100 test



* BRETANOMYCES Selective Agar for wine Dekkera alterative, TPL937, 20 or 200 u



* Campy Selective Preston Agar DMT184, 500g + Suppl.SLM131, for 5 liters

* Raw sterile Matrix (meat, fish, eggs, plants... all you need), for microbiological validations, to inoculate know quantities of different strains, KMT363, 10 x 25g

