


COSMETIKIT®

Laboratorios
MICROKIT



¡Pasión por la creatividad!



COSMETIKIT® is not just a kit for beginners and cover cosmetic microbiology record.

COSMETIKIT® is the whole ready-to-use media that have proved most suitable in cosmetic microbiology, both in its composition as presentation in a format suitable for inclusion in the mass inoculation.

Nobody can make a microbiological control of cosmetics more appropriate and comprehensive than you will with a COSMETIKIT®. Because not follow outdated protocols, but the protocols validated by MICROKIT®, which are at least 140-250% more effective than the Pharmacopeia, ISO standards and derived protocols.



- ★ In counts of aerobes and fungi, you will multiply by 10 the lower limit of quantification, avoiding frequent false-negative results below the guideline value required by the health authorities.
- ★ On detection of the 5 cosmetic pathogens (*Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Candida albicans*, *E.coli-Coliforms* and *Burkholderia cepacia*) you will not risk the dangerous and likely occurrence of false negatives, using the method that ensures maximum sensitivity and specificity.

If you prefer, you can ask us its dehydrated media components:

- * LPT Neutralizing Broth (Lethen updated)
- * Rosa Bengal Caf.Agar
- * Mannitol Salt Agar
- * Mugplus Agar
- * LPT purple Agar
- * Cetrimide Agar
- * Biggy Nickerson Agar
- * BCPT Agar.

For ordering us all of it, ask free our scheme of validated microbiological testing of cosmetics.

- ➔ Also available the new **COSMETIKIT®-EASY**, to save the melting of solid media in the analysis by mass inoculation of your cosmetic products.
- ➔ **NEW! COSMETIKIT®-WATER** for complete microbiological cosmetic use water, without filtering devices and with maximum efficiency.